

Prisma

is a revolutionary thin film micropump

which can be used as an innovative pumping system in wearable insulin delivery devices. Prisma will introduce a breakthrough approach to the treatment of diabetic patients.



Follow Prisma on _____



prisma-horizon.eu

Partners _____




The Prisma Project was funded by the EU Commission in the framework of the Horizon Europe – EIC Transition Open programme. Grant agreement: 101057436

Micro-mechanical pump for
the next generation insulin
delivery systems

Pump devices based on traditional technology

- ✗ Visible
- ✗ Non accurate
- ✗ Complex



Devices based on Prisma micro pumping system


- + Size, way smaller and discrete
- + Higher drug delivery accuracy within the 5% expected range even with the smallest flow
- + Drastic reduction in energy consumption

PAVING THE WAY TOWARDS MULTI HORMONE THERAPIES



How does Prisma thin microfilm pump work?

The electromechanically active materials (thin film) generate a micro-actuation, making the liquid flow forward.



How Prisma will improve insulin therapies?



Discreteness

The reduced size of Prisma allows for a free design of shapes that will result in a pump that the user can forget about



Simplicity

The reduced energy consumption, combined with the accuracy and size of Prisma, will allow for a better closed loop system that will simplify the user experience



Reliability

The precise drug delivery will provide Prisma users with a reliable alternative to insulin pens



Multi-hormone therapy

Thanks to the reduced size of the pump, the multi-hormone therapy will become a reality soon